Gebrauchsanleitung TPP - TubeSpin® Bioreaktor / *Instruction for Use TPP TubeSpin® Bioreactor*

Anleitung

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Instruction

The TubeSpin® Bioreactor is intended for performing media optimization and process development experiments of suspension cells in high numbers, cultivation in shakers at 37 °C temperature and centrifugation.

- *Introduction of medium: put tube into appropriate rack, fill medium as with the classic 50 ml centrifuge tube design using a pipette and allowing the medium to flow down the sides of the tube, please note the recommended filling volume.*
- *Closure: screw down filter screw cap tightly. Holes of partly different size above a PTFE membrane serve as sterile barrier enable optimal gas exchange with the surrounding atmosphere; and minimize liquid loss by evaporation.*
- *Cultivation: mass transfer is achieved by shaking in an appropriate incubation shaker (temperature, CO₂ and humidity controlled; regulation of shaking frequency).*
- *Shaker: Please refer to the relevant shaker information, for ex. www.kuehner.ch*
- *Centrifugation: use the applicable centrifuge rotor or adapter.*
- *Maximum centrifugal acceleration force (gx): can be influenced by following criteria: use of organic solutions, temperature below 0 °C, centrifuging without centrifuge adapter, etc.*

Information to the shanking culture can be found on www.shakingtechnology.com.

Caution

Consult the maximum centrifugal acceleration force (gx or RCF) as well as the necessity to use the applicable centrifuge adapter. Alcohol < 90 % plus rubbing can dissolve and/or blur the white color print on the tube. The TubeSpin® Bioreactors should only be inserted/withdrawn into/from the rack when the shaker is at a completely halt. Tube and cap are not intended to be used in freezing processes.

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Technical Data

Materialien

Schraubkappe	PE
Membrane	PTFE
Röhrchen	PP / PS

Materials

Screw Cap	PE
Membrane	PTFE
Tube	PP

Abmessungen	Measurements	87015	87017	87050	87450	87600
Volumen grad.	Volume grad.	15	15	50	450	600
Länge mm	Length mm	120	120	115	145	181
Durchmesser mm	Diameter mm	17.1	17.1	30	96	98
(g) max	(gx) max	15'500	15'500	15'500	3'500	3'500
Form	Form	Konisch/ conical	Rund/ round	Konisch/ conical	Konisch/ conical	Konisch/ conical
Optimale Einfüllmenge ml	Optimal filling volume ml	1 – 10	1 – 10	1 – 35	>300	>400
Schüttler: Empf.	Shaker:					
Durchmesser mm	Recom. amplitude mm	50	50	50	50	50
Schüttler: Empf.	Shaker:					
Geschwindigkeit rpm	Recom. speed rpm	250	250	180	150	150
Barcode	Barcode	Präfix 7 zur Bestellnummer z.B: 787050 = Bioreaktor 50 mit Barcode Prefix 7 to order number for ex. 757050 = Bioreaktor 50 with barcode				
Schraubkappe "solid"	Screw cap „solid“	--	--	--	87355	87355
Schraubkappe "filter"	Screw cap „filter“	--	--	--	87356	87356