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Instructions

TPP tissue culture flasks are intended for cell culture with necessary sterile handling techniques. The growth area is opto-mechanically treated for optimal cell adhesion and proliferation. The angled neck decreases the risk of media contact with the inside of the screw cap.

- *Open the product packaging and inspect for intact bag and product.*
- *Fill medium: Introduce medium or inoculate cells with a pipette alongside the wall, or direct injection on the growth area of the flask.*
- *Tightly screw closed the filter screw cap with the integrated hydrophobic membrane for sterile gas exchange*
- *Incubate: Lay the flask on the incubator shelves for short or long term storage.*
- **Open flask with peel-off foil:** *Peel-off foil is peeled off before access to the cells. Grasp the flab and peel off the foil. Re-closure is not possible.*
- *Harvest cells with large access from above to flask chamber*
 - *Harvest adherent cells by trypsinization or scrape with the cell scraper.*
 - *Harvest suspension cells by exchange of medium.*

Advice and Recommendation

- *Stable ring enables slip-free stacking of multiple flasks also in combination with the standard TPP tissue culture flasks.*
- *Barrier creates a growth area of 115 cm²*
- *Both flasks are available with Filter screw cap with a hydrophobic PTFE membrane with a pore size of 0.22 µm only. If a «VENT» screw cap is used, the build-up of inside pressure of more than 0.03 bar could press open the flask.*
- *A controlled level base surface and regular use of feeding medium are essential to uniform tissue growth.*
- *Sources of vibrations must be eliminated to avoid circular growth pattern.*

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